

Microcatheter Components' GRIPLOCK™ thin-wall adhesive tubing provides the ultimate flexible assembly solution for layer bonding challenges.

Supplying a bonding layer* for reflowed catheter assemblies, GRIPLOCK™ is offered in an array of standard dimensions with two standard wall thickness options; custom sizes are available by request.

INSIDE DIAMETER		AVERAGE WALL THICKNESS			
		THIN WALL		ULTRA-THIN WALL	
Inches	mm	Inches	mm	Inches	mm
0.010	0.254	0.0015	0.0381	0.001	0.0254
0.015	0.381	0.0015	0.0381	0.001	0.0254
0.020	0.508	0.0015	0.0381	0.001	0.0254
0.025	0.635	0.0015	0.0381	0.001	0.0254
0.030	0.762	0.0015	0.0381	0.001	0.0254
0.035	0.889	0.0015	0.0381	0.001	0.0254
0.040	1.016	0.0015	0.0381	0.001	0.0254
0.045	1.143	0.0015	0.0381	0.001	0.0254
0.050	1.270	0.002	0.0508	0.0015	0.0381
0.055	1.397	0.002	0.0508	0.0015	0.0381
0.060	1.524	0.002	0.0508	0.0015	0.0381
0.065	1.651	0.002	0.0508	0.0015	0.0381
0.070	1.778	0.002	0.0508	0.0015	0.0381
0.075	1.905	0.002	0.0508	0.0015	0.0381
0.080	2.032	0.002	0.0508	0.0015	0.0381
0.085	2.159	0.002	0.0508	0.0015	0.0381
0.090	2.286	0.0025	0.0635	0.002	0.0508
0.095	2.413	0.0025	0.0635	0.002	0.0508
0.100	2.540	0.0025	0.0635	0.002	0.0508
0.105	2.667	0.0025	0.0635	0.002	0.0508
0.110	2.794	0.0025	0.0635	0.002	0.0508
0.115	2.921	0.0025	0.0635	0.002	0.0508
0.120	3.048	0.0025	0.0635	0.002	0.0508
0.125	3.175	0.003	0.0762	0.0025	0.0635
0.130	3.302	0.003	0.0762	0.0025	0.0635
0.135	3.429	0.003	0.0762	0.0025	0.0635
0.140	3.556	0.003	0.0762	0.0025	0.0635
0.145	3.683	0.003	0.0762	0.0025	0.0635
0.150	3.810	0.003	0.0762	0.0025	0.0635

Our customers understand their applications and design requirements best. Suitability for use is their sole responsibility.

©2019 Microcatheter Components